

# OtoClear<sup>®</sup>

Ear Irrigation Tip

A Safer Choice for  
Ear Wax Removal

Made in the USA

For a better patient experience, try the Bionix<sup>®</sup> OtoClear Ear Irrigation Tip.

The OtoClear Tip provides safe ear irrigation with three divergent streams that direct water comfortably to the ear canal wall and away from the ear drum.

Using the OtoClear SprayWash kit costs less than many other irrigation devices and there's no messy backsplash. The OtoClear Tip can be used with these ear irrigation systems from Bionix: the OtoClear SprayWash Kit, OtoClear Aquabot<sup>®</sup> Kit and the Portable Electronic Irrigation System.



**Get your FREE Sample**

See back of card for details.



## Request a FREE OtoClear® Tip sample!

### Four Simple Ways to Request an OtoClear Ear Irrigation Tip sample:

- 1 Scan QR Code or visit **bionix.com/TryOtoClear**  
*International customers – please contact  
your distributor*
- 2 Call: 1.800.551.7096 (+1 419.727.8421)
- 3 Email: [bionix@bionix.com](mailto:bionix@bionix.com)  
[international@bionix.com](mailto:international@bionix.com)



1670 Indian Wood Circle  
Maumee, OH 43537  
United States

800.551.7096 (+1 419.727.8421)  
800.455.5678 (+1 419.727.8426)  
[bionix.com](http://bionix.com)  
[bionix@bionix.com](mailto:bionix@bionix.com)  
[international@bionix.com](mailto:international@bionix.com)



SCAN ME