

Provider Promotion

**Purchase OtoClear®
Ear Irrigation Tips (#7202)
and receive a complimentary
20 pack of OtoClear Tips
(#7202)**

150,000 times a week, someone in the U.S. needs earwax removal. Irrigation is often a patient-preferred method over curettage. Plus, it's effective even for large, impacted wax plugs, and can be easily performed.

Bionix OtoClear Ear Irrigation Tips and Systems are on contract with Vizient and HPG.

Contact Zak Huttner

✉ zhuttner@bionix.com

Part of the **OneLight®**
Illumination Platform



Bionix OtoClear Tips make irrigation an exceptionally safe and gentle procedure, due to unique design features:

3-stream tip directs water to ear canal walls, reducing risk of injury to the tympanic membrane

Large exit portals and ribbing prevent sealing to the ear canal and back-pressure buildup

A flared design prevents over-insertion and fits patients of all ages

Bionix OtoClear Ear Irrigation Tips are the clear choice for safe removal of impacted cerumen.


bionix®
Innovation With Impact